



# The Survivalink® FirstSave® AED



The easy-to-use features of the FirstSave® AED allow virtually anyone with proper training to defibrillate SCA victims quickly and effectively, improving the chances of survival dramatically.

- + Opening the lid automatically turns the FirstSave unit on; it's ready for use.
- + Clear, concise voice-prompts guide the user.
- + RescueReady® feature performs daily, weekly and monthly self-tests.
- + The unit is powered by a five-year or two-year lithium battery, providing maintenance-free operation.
- + Customized defibrillation with escalating, variable STAR™ Biphasic or STAR Monophasic waveforms.

The Rescue data card uses CompactFlash™ technology to record the victim's ECG information for data tracking.

Smartgauge indicator allows user to monitor the battery's status at any time.

Exclusive one-button operation eliminates confusion and uncertainty.

RescueReady indicator turns from green to red, and a loud audio alert sounds if service is required.

Non-polarized, interchangeable electrodes diagram proper placement and reduce confusion at the rescue scene.

Pre-connected electrodes save time and minimize possible errors.



#### Product Code

9200-2 FirstSave AED Biphasic Waveform, complete with Extended Life Lithium Battery and Internal Memory.

## ACCESSORIES AND DATA MANAGEMENT

The Survivalink product line includes a number of accessory products and data management tools to administer an AED program.



### FIRSTSAVE AED TRAINER

This training unit has the same ease-of-use features as the FirstSave AED without the capability to deliver a shock. Using a 9-volt battery, it allows hands-on practice at the student's pace in a variety of simulated rescue scenarios. It also includes an advance button that allows instructors to instantly change the rescue scenario for more realistic training.

#### Product Code

9163 FirstSave AED Training Device



### KIT CARE AED READY KIT

This package provides emergency supplies for initial response with the FirstSave AED. The kit includes items recommended by the AHA, such as a disposable razor, scissors, disposable towel, barrier mask and filter, and disposable latex-free gloves.

#### Product Code

9911 Kit Care AED Ready Kit



### STANDARD ELECTRODES

Survivalink offers single-use, interchangeable electrodes with a large gel surface area to assure maximum success in defibrillation.

#### Product Code

9130 Defibrillation, 2 Year Electrodes



### BATTERIES

Both two- and five-year lithium batteries are available with IntelliSense® technology that monitors the life left in the battery.

#### Product Code

9140 FirstSave Standard 2 Year Lithium Battery  
9141 FirstSave Extended Life 5 Year Lithium Battery



### SOFT-SIDED CARRYING CASE

The soft-sided carrying case is a simple, convenient option for transporting a basic AED ready kit and the FirstSave AED.

#### Product Code

9012 Soft-Sided Carrying/Storage Case

### AED SECURITY WALL CABINET

The AED Security Wall Cabinet is a premium wall-mount storage case that offers additional security features for displaying and storing the FirstSave AED.

#### Product Code

5612 AED Security Wall Cabinet  
24" x 14" x 5" (60.9 cm x 35.6 cm x 12.7 cm)

